

New Product Introduction

TAKE A BREATH™ Portable Handheld Mesh Nebulizer

For Adults and Children

Item # BJ240700

- Super small nebulizer with vibrating mesh technology fits discreetly in the palm of your hand
- Vibrating mesh technology is designed to aerosolize liquid medication into an ultra fine mist for fast and efficient medication delivery
- Medication cup can be removed for easy cleaning
- Nebulizer can be charged with the included USB cable and power adapter - no batteries to replace
- Two nebulization modes: uniform-nebulization mode and micro-wave mode
- Tubeless and cordless design makes nebulizing easy and convenient wherever you go
- Two masks included for use by children and adults
- Small and lightweight - weighs only 3.3 oz and measures 4" long



Take A Breath™ Portable Handheld Mesh Nebulizer

Perfect for Children and Adults

- INCLUDES**
- Portable Mesh Nebulizer
 - Removable Medication Cup
 - Adult Mask
 - Child Mask
 - Mouthpiece
 - USB Charging Cable
 - USB Power Adapter
 - Storage/Travel Carry Case

BLUE JAY

Net. Wt. 1.1 oz (31.9 g)

Item # BJ240700

Retail box



Included accessories:
storage case, adult mask,
child mask, mouthpiece,
USB power adapter & cable

SPECIFICATIONS	
Power Supply	DC 3.7 V (rechargeable battery)
Power Consumption	< 4.0 W
Nebulization Rate	0.15 ml/min~0.9 ml/min
Particle Size	MMAD < 5 µm
Medication Cup Capacity	6 ml
Noise	≤ 50dB
Vibration Frequency	130kHz ± 10%
Weight	3.3 oz (94 g)
Dimensions	2.4" L x 1.6" W x 4" H
Operating Temperature Range	41° F to 104° F
Operating Humidity Range	30% to 85% R.H.
Storage Temperature Range	-4° F to 131° F
Storage Humidity Range	30% to 85% R.H.
Operating Cycle	10 minutes
Battery Life	Up to 70 minutes continuous use
In the Box	Nebulizer, Adult Mask, Child Mask, Mouthpiece, USB Charging Cable, USB Power Adapter, Storage Case
Warranty	Nebulizer - 1 year Mesh Module - 6 months Accessories - 3 months



**Fits discreetly
in the palm of
your hand**