

New Product Introduction



MAKE EATING EASIER™ BIG Grip Utensil Set



Ribbed handles provide enhanced grip

Item	MAKE EATING EASIER™ BIG Grip Utensil Set
Item #	BJ100210
Features	<ul style="list-style-type: none"> • Includes Tablespoon, Fork, Rocker Knife, and FREE Bonus Soup Spoon • Comfortable 1.5" wide built up rubber handle is easy to hold for those with a weak grip • Ribbed handles provide enhanced grip • Stainless steel utensils ensure durability and longevity • Dishwasher safe (up to 125°) • Latex-free handles
Barcode #	615435815619
Case Pack	Case/20



Attractive retail box